

Scavenging Cube

Proper scavenging of waste anesthetic gas can be challenging. Too much vacuum sucks the gas past the patient; not enough causes you to be exposed.

The VetEquip Scavenging Cube addresses these problems like never before.

The Scavenging cube creates low-pressure/high-flow exhaust for chambers such as the VetEquip Vented Chamber, and low-flow/high-pressure vacuum for circuits like the CX-series nosecones or similar. Vent to a VaporGuard filter or in-house evacuation for easy removal of all waste gas.

THE SCAVENGING CUBE ALLOWS USE OF A CIRCUIT AND CHAMBER WITH ANY TYPE OF EVACUATION.

Shown with 2L Vented Chamber and CX-R nosecone (both sold separately)



Scavenging Cube
for Anesthetic Gas

Item Number: 931600



(Shown with optional CX-R and vented induction chamber)

www.kftechnology.it – info@kftechnology.it